



268-S

Befestigungsstück mit Doppelgelenk für Kopfkalotten/Kopfplatten
Connecting fixture with double joint for head rests

Befestigungsstück für Kopflagerungszubehör (Kopfkalotten und Kopfplatten), komplett aus Chromnickelstahl (Edelstahl) gefertigt. Doppelgelenk mit zwei Präzisions-Zahnkranzkupplungen und Selbsthemmung, das über eine zentrale Griffschraube fixiert wird. Griffschraube aus Chromnickelstahl. Gleitschiene im Standardmaß 25 x 10 mm zur Aufnahme verschiedener Kopfkalotten und Kopfplatten. Befestigung an vorhandenen Vierkantaufnahmen 20 x 20 mm am Rückenrahmen des OP-Tisches oder an entsprechenden Verbindungsbügeln mit Vierkantaufnahme zwischen den Gleitschienen des OP-Tisches. Einsteckvierkant des Befestigungsstückes mit automatischer Auszugssperre. Gewicht: 2,5 kg. Stand 01/2010. Technische Änderungen vorbehalten.

Connecting fixture for head positioning accessories (head rests and head plates), completely made of stainless steel. Double joint with 2 crown gear adjustments, self-retarding, to be fixed by a central locking screw. Central locking screw with handle completely made of stainless steel. Side rail with standard dimension 25 x 10 mm for fixing different kinds of head rests and head plates. Connecting fixture to be attached with its adapter square bar 20 x 20 mm to existing square tubes of the operating table's upper back plate frame or to head rest attachments bars with square tubes connecting the side rails of the operating table's upper back plate. Adapter bar with automatic extension limiter. Weight: 2,5 kg. Release 01/2010. Subject to technical alterations.

Schindler
Krankenhausentwicklung
GmbH

Postfach 15 08, 22905 Ahrensburg
Kurt-Fischer-Str. 27 c, 22926 Ahrensburg
Deutschland/Germany
T. +49 (0) 41 02 / 80 6 80 - 0
F. +49 (0) 41 02 / 80 6 80 - 28
E. info@schindler-krankenhausentwicklung.de
www.schindler-krankenhausentwicklung.de